

Figure 1. M-mode and speckle tracking overview. (A) M-mode: Lower panel represents tissue displacement passing through the M-mode beam over time. Vertical white line indicates diaphragm displacement. (B) M-mode does not track the same region of tissue over time due to simultaneous lateral movement. The solid curved line represents initial diaphragm position. During breathing, there is cephalocaudad excursion of tissue represented by dashed curved lines and lateral diaphragm movement illustrated by the geometric shapes shifting from right to left. The M-mode beam (large gray arrow) passes through different points on the tissue during breathing, thus, displacement measurement (small black arrow) represents different points on the tissue. (C) Speckle tracking: B-mode ultrasonography in the left panel shows regions of interest (ROI) on the diaphragm tissue (small white boxes, indicated by white arrowheads). Movement of the ROIs in the X and Y planes are shown in the right upper and lower panels, respectively. See Supplemental Figure S1 for additional details. (D) Speckletracking measures displacement in 2 dimensions and therefore can better characterize true motion. The actual measurement (black arrows) is fixed relative to tissue position.

152x115mm (300 x 300 DPI)



This is the author manuscript accepted for publication and has undergone full peer review but has not been through the copyediting, typesetting, pagination and proofreading process, which may lead to differences between this version and the Version record. Please cite this article as doi:10.1002/mus.25380.

Title: Speckle tracking as a method to measure hemidiaphragm excursion

Counts: Paper: 1080 words; Abstract: 146 words; Title: 63 characters; Tables: 0 (1 supplemental);

Figures: 1 (1 supplemental); Videos 4

Authors: Stephen A. Goutman, MD, MS¹, James D. Hamilton, PhD², Blake Swihart, MA¹, Bradley Foerster, MD, PhD^{3,4}, Eva L. Feldman, MD, PhD^{1,5}, Jonathan M. Rubin, MD, PhD³

Affiliations:

¹Department of Neurology

²Department of Internal Medicine

³Department of Radiology

⁴Ann Arbor Veterans Administration Healthcare System, Ann Arbor, MI

⁵A. Alfred Taubman Medical Institute University of Michigan, Ann Arbor, MI

Study Funding: A. Alfred Taubman Medical Institute

Corresponding Author:

Stephen A Goutman, MD, MS Department of Neurology F2647 UH SOUTH, SPC 5223 1500 E Medical Center Dr Ann Arbor, MI 48109-5223 Phone: 734-936-8586

Fax: 734-936-5185

ACCE

e-mail: sgoutman@med.umich.edu

Running Title: Diaphragm speckle tracking



Abstract:

Introduction: Diaphragm excursion measured via ultrasound may be an important imaging outcome measure of respiratory function. We developed a new method of measuring diaphragm movement and compared this to the more traditional M-mode method.

Methods: Ultrasound images of the right and left hemidiaphragms were collected to compare speckle tracking and M-mode measurements of diaphragm excursion. Speckle tracking was performed in EchoInsight (Epsilon Imaging, Ann Arbor, MI).

Results: Six healthy subjects without a history of pulmonary diseases composed this proof of concept study. Speckle tracking of the diaphragm is technically possible. Unlike M-mode, speckle tracking carries the advantage of reliable visualization and measurement of the left hemidiaphragm.

Discussion: Speckle tracking accounted for diaphragm movement simultaneously in the cephalocaudad and mediolateral directions, unlike M-mode which is one-dimensional. Diaphragm speckle tracking may represent a new and more robust method for measuring diaphragm excursion, especially for the left hemidiaphragm.

Key Words:

- 1. Diaphragm ultrasound
- 2. Diaphragm excursion
- 3. Imaging biomarker
- 4. Speckle tracking
- 5. M-mode
- 6. Respiratory imaging biomarker

Accepted Ar

Introduction

M-mode ultrasonography can quantify diaphragm excursion,¹⁻⁴ although there are limitations: (1) the left hemidiaphragm is rarely visible due to the smaller acoustic window in the left upper quadrant;³⁻⁵ and (2) M-mode is a one-dimensional (1D) measurement of tissue passing through the ultrasound beam. Therefore, although tissue moving towards and away from the ultrasound probe can be measured, if moves laterally concurrently, this vector cannot be captured and consequently the M-mode does not reflect true motion.

Speckle tracking is used in echocardiography^{6, 7} and can evaluate organ motion.⁸ Speckle is created when ultrasound signals are reflected from microstructure within tissue and other scattering media. Because these micro-scatterers are in slightly different positions in an ultrasound beam, their contributions to the backscattered ultrasound signal may be positive or negative depending on the phase of the ultrasound wave at each scatterer position. The result is an overlaying pattern of light and dark (referred to as salt and pepper) textured appearance in ultrasound images. This pattern, which looks like noise, is often stable within small displacements, and can be tracked during motion using computer-based image processing techniques. Such tracking can be performed in 2- or 3-dimensions (2D or 3D), in contrast to M-mode, which instead measures 1D tissue motion and does not track a fixed region of tissue. We hypothesized that speckle tracking robustly measures diaphragm movement, as it accounts for the 2D motion vector rather than 1D motion with M-mode imaging.

Materials and Methods

Participants. Healthy participants had no history of pulmonary or neuromuscular disease. All subjects provided written informed consent, and the study received Institutional Review Board approval. Age and gender were recorded.

Ultrasound protocol. Ultrasounds were performed in a semirecumbent position using a Philips Epic ultrasound machine (Philips Ultrasound, Bothel, WA) with a standard abdominal imaging curved array transducer, C5-1. For speckle tracking loops, the probe was oriented coronally at the mid-axillary line, and each corresponding hemidiaphragm was imaged between the ribs. Ultrasound images recorded the right and left hemidiaphragm during both normal and deep breathing. For M-modes, the same probe was oriented sagittally subcostally. The M-mode sampling line was directed so that the motion would be oriented directly toward the probe, i.e. in a cephalocaudad direction. At least 3 B- and M-mode samples for the right and left hemidiaphragms during both normal and deep inspiration were performed. M-mode displacements were averaged for analysis.

Speckle Tracking. Speckle tracking was performed using EchoInsight (Epsilon Imaging, Ann Arbor, MI), a commercial software package designed to track 2D speckle motion in cardiac echograms and other applications. Following import of DICOM files into EchoInsight, a region of interest (ROI) curve was placed over the diaphragm, utilizing the *define wall* tool (Figure 1, Videos S1-S4, Supplemental Figure). The ROI width was 0.5cm. As shown in the figure, EchoInisight creates 6 individual measurement regions, or segments, from the ROI curve. To create evenly spaced ROI segments, the center of the ROI curve was set visually using control+left mouse click. The image loop was reviewed to ensure the ROI did not drift from the diaphragm tissue, and if so, the tracking operating mode was switched (see below). Each ROI curve measured cephalocaudad displacement, defined as the square root of the X displacement squared plus the Y displacement squared, averaged over the 6 ROI segments for each

loop. For analysis, the diaphragm was visualized for a complete respiratory cycle; if a shadow obscured the speckle over the diaphragm, the image was excluded from analysis.

The Echolnsight speckle tracking algorithm was modified to have 2 operating modes, which controlled the resolution of the image data prior to tracking. In fine resolution mode, B-mode images were spatially filtered using a 2D, 5 point Hanning filter prior to tracking. Coarse resolution mode used a 15 point Hanning filter. The 2 modes allowed the operator of Echolnsight to trade off between robustness and precision of the tracking algorithm by adjusting the spatial frequency content of the tracked data. The fine resolution algorithm was used on diaphragm loops with a sharp and clearly visualized diaphragm border. For noisier data and image loops with a diffuse appearance of the diaphragm border, the coarse tracking algorithm was used. The operator determined which mode to implement by visually assessing the tracking of the ROI curve throughout the respiratory cycle.

Statistical analysis. All analyses were performed with SAS 9.4 (SAS Institute, Carey, NC).

Results

Six subjects (3 women) underwent ultrasound assessment. Mean age was 57.7 years with a standard deviation of 9.1 years. M-mode ultrasonography and speckle tracking each visualized all 6 subjects' right hemidiaphragms during both normal and deep inspiration. On the left, M mode captured 2 subjects' hemidiaphragms during normal inspiration and none during deep inspiration, whereas speckle tracking captured all 6 hemidiaphragms during normal breathing and 5 of 6 during deep breathing.

Supplemental Table S1 (available online) highlights displacements by subject. The average difference between right speckle tracking and M-mode for normal and deep inspiration was 0.30 cm (n=6) and -

0.65 cm (n=6) respectively and on the left was 1.00 cm (n=2) for normal inspiration. No M-modes were available on the left for comparison during deep inspiration.



We provide an initial description of speckle tracking for evaluating diaphragm motion. This technique has 2 distinct advantages. First, it more accurately evaluates diaphragm movement compared to Mmode, because speckle tracking measures the motion of a defined point or region on the diaphragm through the image loop, whereas M-mode records motion of variable regions of tissue passing through the M-mode line. Second, speckle tracking can quantify left hemidiaphragm displacement, a nearly impossible task with M-mode.

One limitation of speckle tracking is additional time needed for software processing of ultrasound studies. In some circumstances, displacement magnitudes were larger for M-mode compared to speckle tracking. The diaphragm does not move uniformly, and speckle tracking can measure a larger region of tissue, including segments without peak displacements. In contrast, M-mode focuses on peak displacement. Analysis of individual speckle segments showed displacement variability of up to 0.5 cm, suggesting that the speckle tracking averages out motion over a larger range of tissue. Speckle tracking can also capture peak displacements; future studies should focus on whether bulk or peak motion is more clinically meaningful. Importantly, as speckle tracking clearly demonstrates displacements in the X and Y directions, we show important limitations of M-mode for diaphragm displacement measurements. If speckle tracking is further validated, it may prove to be a more useful and reliable method for longitudinal assessment of diaphragm movement with ultrasound compared to current techniques.

Acknowledgements

The authors thank Danielle Dobbs and Vanessa Allen for their assistance with medical illustrations.

Research funding support was provided by the A. Alfred Taubman Institute. The authors also thank

Epsilon Imaging for providing software and technical support. DICOM image transfer and anonymization

done by Clinical Trial Processor (Radiological Society of North America, Chicago, IL).

Abbreviations

1D: one dimensional

2D: two dimensional

3D: three dimensional

ROI: regional of interest

References

- 1. Sarwal A, Walker FO, Cartwright MS. Neuromuscular ultrasound for evaluation of the diaphragm. Muscle & nerve 2013;47:319-329.
- 2. McCool FD, Tzelepis GE. Dysfunction of the Diaphragm. New England Journal of Medicine 2012;366:932-942.
- 3. Gerscovich EO, Cronan M, McGahan JP, Jain K, Jones CD, McDonald C. Ultrasonographic evaluation of diaphragmatic motion. Journal of ultrasound in medicine: official journal of the American Institute of Ultrasound in Medicine 2001;20:597-604.
- 4. Boussuges A, Gole Y, Blanc P. Diaphragmatic motion studied by m-mode ultrasonography: methods, reproducibility, and normal values. Chest 2009;135:391-400.
- 5. Scott S, Fuld JP, Carter R, McEntegart M, MacFarlane NG. Diaphragm ultrasonography as an alternative to whole-body plethysmography in pulmonary function testing. Journal of ultrasound in medicine: official journal of the American Institute of Ultrasound in Medicine 2006;25:225-232.
- 6. Blessberger H, Binder T. NON-invasive imaging: Two dimensional speckle tracking echocardiography: basic principles. Heart 2010;96:716-722.
- 7. Blessberger H, Binder T. Two dimensional speckle tracking echocardiography: clinical applications. Heart 2010;96:2032-2040.
- 8. Rubin JM, Feng M, Hadley SW, Fowlkes JB, Hamilton JD. Potential use of ultrasound speckle tracking for motion management during radiotherapy: preliminary report. Journal of ultrasound in medicine: official journal of the American Institute of Ultrasound in Medicine 2012;31:469-481.

Figures

Figure 1. M-mode and speckle tracking overview. (A) M-mode: Lower panel represents tissue displacement passing through the M-mode beam over time. Vertical white line indicates diaphragm displacement. (B) M-mode does not track the same region of tissue over time due to simultaneous lateral movement. The solid curved line represents initial diaphragm position. During breathing, there is cephalocaudad excursion of tissue represented by dashed curved lines and lateral diaphragm movement illustrated by the geometric shapes shifting from right to left. The M-mode beam (large gray arrow) passes through different points on the tissue during breathing, thus, displacement measurement (small black arrow) represents different points on the tissue. (C) Speckle tracking: B-mode ultrasonography in the left panel shows regions of interest (ROI) on the diaphragm tissue (small white boxes, indicated by white arrowheads). Movement of the ROIs in the X and Y planes are shown in the right upper and lower panels, respectively. See Supplemental Figure S1 (available online) for additional details. (D) Speckle tracking measures displacement in 2 dimensions and therefore can better characterize true motion. The actual measurement (black arrows) is fixed relative to tissue position.

