

Table S2. New primers developed for amplification and sequencing of the complete cytochrome B and for internal sequencing of ND2. Previously published primers (Kocher 1995) were also used for amplification and sequencing of ND2.

Gene	Primer	Direction	Primer sequence 5'-3'	Location
Cyt B	14439.chls	F	TTGAAAAACCACCGTTGTTATTCAACTA	Flanking
	tR.chls	R	CCATTTTCCGGTTTACAAGACCG	Flanking
	523.cytBint	F	CTAACCCGCTTCTTCGCCTTC	Internal
	238.cytBint	R	GCCGTTAGCGTCCATATTTTCG	Internal
	637.cytBint	R	GAGATTTTGTTCGGCGTCTGAG	Internal
ND2	412.ND2int	F	CCCTTCGCCCTAATTCTTCAAC	Internal
	597.ND2int	R	GAGGTAGGTTAGAAGGGTAAGGAGGGTG	Internal